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
Insertion image du pied

GESPODO

The screenshot shows the main window of the 'GesPodo Design' application. On the left, a vertical toolbar contains icons for file operations, editing, and production. A 'Vue' (View) panel is open, showing options for displaying the foot model and its components. On the right, a vertical toolbar contains icons for various therapeutic modifications. At the bottom right, there are icons for undo/history and image insertion. The 'GESPODO' logo is visible in the top right and bottom right corners.

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

Menu torsion

The screenshot displays the GESPODO software interface. On the left, a vertical toolbar contains various icons, with the 'Torsion' icon (a foot with a curved arrow) circled in blue. A blue arrow points from this icon to a text box. The central area shows a 3D model of a foot in green, positioned on a grid. Two blue lines are drawn on the sole of the foot, indicating the torsion axis. On the right, the 'Zone de torsion' settings panel is open, showing options for 'type' (Torsion globale, Torsion d'avant-pied, Torsion d'arrière-pied), 'Coté' (Interne, Latéral), and numerical values for 'Angle (°)', 'Position de l'axe vertical (mm)', 'Localisation de l'axe', 'Position de la semelle / sol', 'Déplacer vers le haut / bas (mm)', and 'Basculer vers l'avant / vers l'arrière (0,0,0)'. The 'OK' and 'Annuler' buttons are at the bottom of the panel. The GESPODO logo is visible in the bottom right corner of the software window.

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
The screenshot displays the GESPODO software interface. On the left, a 3D model of a pink shoe sole is shown with various dimensions: 285, 153, 62, 76, 90, 186, h: 31, h: 32, z: 7, z: 4,4, and z: 3,6. A 'Modifications' dialog box is open in the center, showing a table of changes. The table has columns for 'ON', 'Type', and 'Description'. The 'ON' column contains checkboxes, the 'Type' column contains text, and the 'Description' column contains text. The 'Modifications' dialog box has buttons for 'Editer', 'Retirer', and 'Rafraichir'. The GESPODO logo is visible in the top right and bottom right corners of the interface. A blue arrow points to a circular icon in the bottom right toolbar, which is highlighted by a blue circle.

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In the center, a 3D model of a foot is shown with various measurement points and dimensions: 285, 90, 153, 62, 186, h: 31, h: 32, z: 4,4, and z: 3,6.

On the right, a 'Paramètres des images' (Image Parameters) window is open, showing a list of image options. The 'Pied amputé' option is selected. Below the list, there is a preview of a foot image and a 'Chercher...' (Search) button.

Text annotations on the right side of the interface provide keyboard shortcuts:

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